

## Fragment RNA

- 1 Dilute the total RNA in nuclease-free ultrapure water to 8.5 µl in the DFP plate.
- 2 Add 8.5 µl EPH.
- 3 Shake at 1600 rpm for 20 seconds.
- 4 Centrifuge at 280 × g for 1 minute.
- 5 Place on the thermal cycler and run the Elution 2-Frag-Prime program.
- 6 Remove from thermal cycler when it reaches 4°C and centrifuge briefly.

## Synthesize First Strand cDNA

- 1 Add 50 µl SuperScript II to FSA. Pipette to mix, and then centrifuge briefly.
- 2 Add 8 µl SuperScript II and FSA mixture.
- 3 Shake at 1600 rpm for 20 seconds.
- 4 Place on the thermal cycler and run the Synthesize 1st Strand program.
- 5 When the thermal cycler reaches 4°C, remove the DFP plate from the thermal cycler.

## Synthesize Second Strand cDNA

- 1 Add 5 µl RSB.
- 2 Add 20 µl SMM.
- 3 Shake at 1600 rpm for 20 seconds.
- 4 Centrifuge at 280 × g for 1 minute.
- 5 Place on the thermal cycler and incubate at 16°C for 1 hour.
- 6 Place on the bench and let stand to bring to room temperature.
- 7 Add 90 µl AMPure XP beads to the CCP plate.
- 8 Transfer all to the CCP plate.
- 9 Shake at 1800 rpm for 2 minutes.
- 10 Incubate at room temperature for 5 minutes.
- 11 Centrifuge at 280 × g for 1 minute.
- 12 Place on a magnetic stand until liquid is clear.
- 13 Remove and discard 135 µl supernatant.
- 14 Leave the CCP plate on the magnetic stand when performing the following wash step.
- 15 Wash 2 times with 200 µl 80% EtOH.
- 16 Use a 20 µl pipette to remove residual EtOH.
- 17 Air-dry for 5 minutes.
- 18 Remove from the magnetic stand.
- 19 Add 17.5 µl RSB.
- 20 Shake at 1800 rpm for 2 minutes.
- 21 Incubate at room temperature for 2 minutes.
- 22 Centrifuge at 280 × g for 1 minute.
- 23 Place on a magnetic stand until liquid is clear.
- 24 Transfer 15 µl supernatant to the ALP plate.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to seven days.

### Adenylate 3' Ends

- 1 Add 2.5 µl RSB.
- 2 Add 12.5 µl ATL.
- 3 Shake at 1800 rpm for 2 minutes.
- 4 Centrifuge at 280 × g for 1 minute.
- 5 Place on the 37°C microheating system for 30 minutes.
- 6 Move to the 70°C microheating system for 5 minutes.
- 7 Place on ice for 1 minute.

### Ligate Adapters

- 1 Add the following.
  - ▶ RSB (2.5 µl)
  - ▶ LIG (2.5 µl)
  - ▶ RNA adapters (2.5 µl)
- 2 Shake at 1800 rpm for 2 minutes.
- 3 Centrifuge at 280 × g for 1 minute.
- 4 Place on the 30°C microheating system for 10 minutes, and then place on ice.
- 5 Add 5 µl STL.
- 6 Shake at 1800 rpm for 2 minutes.
- 7 Centrifuge at 280 × g for 1 minute.
- 8 Perform steps 9 through 24 using the **Round 1** volumes.
- 9 Add AMPure XP beads.

|                 | Round 1 | Round 2 |
|-----------------|---------|---------|
| AMPure XP beads | 42 µl   | 50 µl   |

- 10 Shake at 1800 rpm for 2 minutes.
- 11 Incubate at room temperature for 5 minutes.
- 12 Centrifuge at 280 × g for 1 minute.
- 13 Place on a magnetic stand until liquid is clear.
- 14 Remove and discard all supernatant.
- 15 Wash 2 times with 200 µl 80% EtOH.
- 16 Use a 20 µl pipette to remove residual EtOH.
- 17 Air-dry for 5 minutes.
- 18 Remove from the magnetic stand.
- 19 Add RSB.

|     | Round 1 | Round 2 |
|-----|---------|---------|
| RSB | 52.5 µl | 22.5 µl |

- 20 Shake at 1800 rpm for 2 minutes.
- 21 Incubate at room temperature for 2 minutes.
- 22 Centrifuge at 280 × g for 1 minute.
- 23 Place on a magnetic stand until liquid is clear.

- 24 Transfer 50 µl supernatant to the CAP plate.
- 25 Repeat steps 9 through 24 with the new plate using the **Round 2** volumes.
- 26 Transfer 20 µl supernatant to the PCR plate.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to seven days.

## Enrich DNA Fragments

- 1 Place on ice and add 5 µl PPC.
- 2 Add 25 µl PMM.
- 3 Shake at 1600 rpm for 20 seconds.
- 4 Centrifuge at 280 × g for 1 minute.
- 5 Place on the thermal cycler and run the PCR program.
- 6 Add 50 µl AMPure XP beads.
- 7 Transfer contents from PCR plate to the CPP plate.
- 8 Shake at 1800 rpm for 2 minutes.
- 9 Incubate at room temperature for 5 minutes.
- 10 Centrifuge at 280 × g for 1 minute.
- 11 Place on a magnetic stand until liquid is clear.
- 12 Remove and discard all supernatant.
- 13 Wash 2 times with 200 µl 80% EtOH.
- 14 Use a 20 µl pipette to remove residual EtOH.
- 15 Air-dry for 5 minutes.
- 16 Remove from the magnetic stand.
- 17 Add 17.5 µl RSB.
- 18 Shake at 1800 rpm for 2 minutes.
- 19 Incubate at room temperature for 2 minutes.
- 20 Centrifuge at 280 × g for 1 minute.
- 21 Place on a magnetic stand until liquid is clear.
- 22 Transfer 15 µl supernatant to the TSP1 plate.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to seven days.

## Check Libraries

- 1 If using a Standard Sensitivity NGS Fragment Analysis Kit:
  - a Dilute the DNA library 1:1 with RSB.
  - b Run 1 µl diluted DNA library.
- 2 If using a DNA 1000 chip, run 1 µl undiluted DNA library.
- 3 Check the size and purity of the sample. Expect the final product to be a band at ~260 bp.

## First Hybridization of Probes

- 1 Add the following to the REH1 plate.
  - ▶ DNA library sample or pool (45 µl)
  - ▶ CT3 (50 µl)
  - ▶ CEX (5 µl)
- 2 For single-plex, add the following to the REH1 plate:
  - ▶ DNA library sample 200ng (11.25 µl)
  - ▶ CT3 (12.5 µl)
  - ▶ CEX (1.25 µl)
- 3 Shake at 1200 rpm for 1 minute.
- 4 Centrifuge at 280 × g for 1 minute.
- 5 Place on the thermal cycler and run the RNA HYB program.
  - ▶ For 2-4 plex pools, each well contains 100 µl.
  - ▶ For single-plex pools, each well contains 25 µl.
- 6 Remove from the thermal cycler immediately after the 90-minute incubation.

## First Capture of Hybridized Probes

- 1 Centrifuge at 280 × g for 1 minute.
- 2 Transfer all to the REW1 plate.
- 3 Add 250 µl (62.5 µl for 1-plex) SMB.
- 4 Shake at 1200 rpm for 5 minutes.
- 5 Incubate at room temperature for 25 minutes.
- 6 Centrifuge at 280 × g for 1 minute.
- 7 Place on a magnetic stand until liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Remove from the magnetic stand.
- 10 Wash two times with 200 µl EWS.
- 11 Mix 28.5 µl EE1 and 1.5 µl HP3, and then vortex to mix.
- 12 For 1-plex, create elution premix in a 1.7 microcentrifuge tube, and then vortex to mix.
  - ▶ EE1 (9.5 µl)
  - ▶ HP3 (0.5 µl)
- 13 Add 23 µl (10 µl for 1-plex) elution premix.
- 14 Shake at 1800 rpm for 2 minutes.
- 15 Incubate at room temperature for 2 minutes.
- 16 Centrifuge at 280 × g for 1 minute.
- 17 Place on a magnetic stand until liquid is clear.
- 18 Transfer 21 µl (9 µl for 1-plex) supernatant to the REH2 plate.
- 19 Add 4 µl (1.7 µl for 1-plex) ET2.
- 20 Shake at 1200 rpm for 1 minute.
- 21 Centrifuge at 280 × g for 1 minute.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to seven days.

## Second Hybridization of Probes

- 1 Add the following.
  - ▶ RSB (20 µl)
  - ▶ CT3 (50 µl)
  - ▶ CEX (5 µl)
- 2 For 1-plex, add the following.
  - ▶ RSB (0.55 µl)
  - ▶ CTE (12.5 µl)
  - ▶ CEX (1.25 µl)
- 3 Shake at 1200 rpm for 1 minute.
- 4 Centrifuge at 280 × g for 1 minute.
- 5 Place on the thermal cycler and run the RNA HYB program.
- 6 Remove from the thermal cycler immediately after the 90-minute incubation.

## Second Capture of Hybridized Probes

- 1 Centrifuge at 280 × g for 1 minute.
- 2 Transfer all supernatant to the REW2 plate.
- 3 Add 250 µl (62.5 µl for 1-plex) SMB.
- 4 Shake at 1200 rpm for 5 minutes.
- 5 Incubate at room temperature for 25 minutes.
- 6 Centrifuge at 280 × g for 1 minute.
- 7 Place on a magnetic stand until liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Remove from the magnetic stand.
- 10 Wash two times with 200 µl EWS.
- 11 Mix 28.5 µl EE1 and 1.5 µl HP3, and then vortex to mix.
- 12 For 1-plex, create elution premix in a 1.7 ml microcentrifuge tube, and then vortex to mix.
  - ▶ EE1 (9.5 µl)
  - ▶ HP3 (0.5 µl)
- 13 Add 23 µl (10 µl for 1-plex) elution premix.
- 14 Shake at 1800 rpm for 2 minutes.
- 15 Incubate at room temperature for 2 minutes.
- 16 Centrifuge at 280 × g for 1 minute.
- 17 Place on a magnetic stand until liquid is clear.
- 18 Transfer 21 µl (9 µl for 1-plex) supernatant to the REC1 plate.
- 19 Add 4 µl (1 µl for 1-plex) ET2.
- 20 Shake at 1800 rpm for 1 minute.
- 21 Centrifuge at 280 × g for 1 minute.

## Clean Up Captured Library

- 1 Add 45 µl (18 µl for 1-plex) AMPure XP beads.
- 2 Shake at 1800 rpm for 1 minute.
- 3 Incubate at room temperature for 5 minutes.
- 4 Centrifuge at 280 × g for 1 minute.
- 5 Place on a magnetic stand until liquid is clear.
- 6 Remove and discard all supernatant.
- 7 Wash 2 times with 200 µl 80% EtOH.
- 8 Use a 20 µl pipette to remove residual EtOH.
- 9 Air-dry for 5 minutes.
- 10 Remove from the magnetic stand.
- 11 Add 27.5 µl RSB.
- 12 Shake at 1800 rpm for 1 minute.
- 13 Incubate at room temperature for 2 minutes.
- 14 Centrifuge at 280 × g for 1 minute.
- 15 Place on a magnetic stand until liquid is clear.
- 16 Transfer 25 µl supernatant to the REA plate.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to seven days.

## Amplify Enriched Library

- 1 Remove seal and add 5 µl PPC.
- 2 Add 20 µl EPM.
- 3 Shake at 1200 rpm for 1 minute.
- 4 Centrifuge at 280 × g for 1 minute.
- 5 Place the midi plate on the thermal cycler and run the EPM AMP program.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days.

## Clean Up Amplified Enriched Library

- 1 Centrifuge at 280 × g for 1 minute.
- 2 Transfer 50 µl to the REC2 plate.
- 3 Add 90 µl AMPure XP beads.
- 4 Shake RAC2 at 1800 rpm for 1 minute.
- 5 Incubate at room temperature for 5 minutes.
- 6 Centrifuge at 280 × g for 1 minute.
- 7 Place on a magnetic stand until liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Wash 2 times with 200 µl 80% EtOH.
- 10 Use a 20 µl pipette to remove residual EtOH.
- 11 Air-dry on the magnetic stand for 5 minutes.
- 12 Remove from the magnetic stand.
- 13 Add 32 µl (8 µl for 1-plex) RSB.
- 14 Shake at 1800 rpm for 1 minute.
- 15 Incubate at room temperature for 2 minutes.
- 16 Centrifuge at 280 × g for 1 minute.
- 17 Place on a magnetic stand until liquid is clear.
- 18 Transfer 30 µl (7.5 µl for 1-plex) supernatant to the REL plate.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to seven days.

## Check Enriched Libraries

- 1 Quantify the libraries according to the [Illumina Sequencing Library qPCR Quantification Guide \(document # 11322363\)](#) or fluorometric method.
- 2 Load 1 µl of the post-enriched library on one of the following:
  - ▶ Advanced Analytical Technologies Standard Sensitivity NGS Fragment Analysis Kit
  - ▶ Agilent High Sensitivity DNA Chip
- 3 Check the size of the library for a distribution of DNA fragments with a size range from ~200 bp–1 kb.
- 4 Proceed to cluster generation.

## Acronyms

| Acronym | Definition                       |
|---------|----------------------------------|
| ALP     | Adapter Ligation Plate           |
| ATL     | A-Tailing Mix                    |
| CAP     | Clean Up ALP Plate               |
| CCP     | cDNA Clean Up Plate              |
| CEX     | Coding Exome Oligos              |
| CPP     | Clean Up PCR Plate               |
| CT3     | Capture Target Buffer 3          |
| DFP     | Depleted RNA Fragmentation Plate |
| EE1     | Enrichment Elution Buffer 1      |
| EPH     | Elute, Prime, Fragment High Mix  |
| ET2     | Elute Target Buffer 2            |
| EWS     | Enrichment Wash Solution         |
| FSA     | First Strand Synthesis Act D Mix |
| HP3     | 2N NaOH                          |
| IEM     | Illumina Experiment Manager      |
| LIG     | Ligation Mix                     |
| LRM     | Local Run Manager                |
| PCR     | Polymerase Chain Reaction Plate  |
| PMM     | PCR Master Mix                   |
| PPC     | PCR Primer Cocktail              |
| REA     | RNA Exome Amplification Plate    |
| REC1    | RNA Exome Clean Up Plate 1       |
| REC2    | RNA Exome Clean Up Plate 2       |
| REH1    | RNA Exome Hyb Plate 1            |
| REH2    | RNA Exome Hyb Plate 2            |
| REL     | RNA Exome Library Plate          |
| REW1    | RNA Exome Wash Plate 1           |

| Acronym | Definition                       |
|---------|----------------------------------|
| REW2    | RNA Exome Wash Plate 2           |
| RSB     | Resuspension Buffer              |
| SMB     | Streptavidin Magnetic Beads      |
| SMM     | Second Strand Marking Master Mix |
| STL     | Stop Ligation Buffer             |
| TSP     | Target Sample Plate              |