

Local Run Manager

PCR Amplicon Analysis Module v1.0

Release Notes

Introduction

These Release Notes detail notable items for the Local Run Manager PCR Amplicon Analysis Module v1.0 release.

Please note that the PCR Amplicon Analysis Module v1.0 requires Local Run Manager Framework v2.0.0 or higher and is not compatible with Local Run Manager Framework v1.3.1 or lower.

For more information about this analysis module and how to use it, refer to the Local Run Manager PCR Amplicon Analysis Module Workflow Guide, available from the Illumina Local Run Manager Support Page, Documentation & Literature.

http://support.illumina.com/sequencing/sequencing_software/local-run-manager/documentation.html

NEW FEATURES:

- Initial Release

DEFECT REPAIRS:

- N/A

KNOWN ISSUES:

- Workflow does not support human hg38 genome due to a known issue related to the configuration of the hg38 build of the human genome.
- When the user selects certain library prep kits, like Nextera DNA, the I7 and I5 index fields will not be marked as required fields.
- Custom Primers don't update in the Local Run Manager run if they are updated in the Control Software.
- If the Local Run Manager Analysis Service is restarted, this could cause a running analysis to crash.

OTHER:

- When using a custom kit, you must always enter your i5 indexes in the MiSeq or forward orientation. Local Run Manager will automatically reverse complement the indexes when writing the sample sheet used for analysis.