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シングルセルシーケンス : [jp.illumina.com/techniques/sequencing/rna-sequencing/ultra-low-input-single-cell-rna-seq.html](https://jp.illumina.com/techniques/sequencing/rna-sequencing/ultra-low-input-single-cell-rna-seq.html)

NovaSeq 6000システム : [jp.illumina.com/systems/sequencing-platforms/novaseq.html](https://jp.illumina.com/systems/sequencing-platforms/novaseq.html)

NextSeq 1000およびNextSeq 2000システム : [jp.illumina.com/systems/sequencing-platforms/nextseq-1000-2000.html](https://jp.illumina.com/systems/sequencing-platforms/nextseq-1000-2000.html)

NextSeq 550システム : [jp.illumina.com/systems/sequencing-platforms/nextseq.html](https://jp.illumina.com/systems/sequencing-platforms/nextseq.html)

Chromium Single Cell Gene Expression : [10xgenomics.com/products/single-cell-gene-expression](https://10xgenomics.com/products/single-cell-gene-expression)

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1. Park JE, Botting RA, Domínguez Conde C, et al. [A cell atlas of human thymic development defines T cell repertoire formation.](#) *Science*. 2020;367(6480):eaay3224. doi:10.1126/science.aay3224
2. Wang R, Dang M, Harada K, et al. [Single-cell dissection of intratumoral heterogeneity and lineage diversity in metastatic gastric adenocarcinoma.](#) *Nat Med*. 2021;27(1):141-151. doi:10.1038/s41591-020-1125-8
3. 10x Genomics. Single cell gene expression support. [support.10xgenomics.com/single-cell-gene-expression](https://support.10xgenomics.com/single-cell-gene-expression). Accessed January 11, 2022.
4. Li YR, Lai HW, Huang HH, Chen HC, Fugmann SD, Yang SY. [Trajectory mapping of the early \*Drosophila\* germline reveals controls of zygotic activation and sex differentiation.](#) *Genome Res*. 2021;31(6):1011-1023. doi:10.1101/gr.271148.120
5. Pelligia J, Münch D, Meneses-Giles P, et al. [Adaptive cell invasion maintains lateral line organ homeostasis in response to environmental changes.](#) *Dev Cell*. 2021;56(9):1296-1312.e7. doi:10.1016/j.devcel.2021.03.027
6. 10x Genomics. Chart your own course with a low-throughput kit for pilot & optimization studies. [pages.10xgenomics.com/rs/446-PBO-704/images/10x\\_LIT000120\\_Product\\_Sheet\\_LT\\_Kit\\_Letter\\_Digital.pdf](https://pages.10xgenomics.com/rs/446-PBO-704/images/10x_LIT000120_Product_Sheet_LT_Kit_Letter_Digital.pdf). Accessed January 3, 2022.
7. 10x Genomics. Explore cellular diversity, cell by cell. [pages.10xgenomics.com/rs/446-PBO-704/images/10x\\_LIT093\\_Product-Sheet\\_Explore-Cellular-Diversity\\_Letter\\_digital.pdf](https://pages.10xgenomics.com/rs/446-PBO-704/images/10x_LIT093_Product-Sheet_Explore-Cellular-Diversity_Letter_digital.pdf). Accessed January 3, 2022.
8. 10x Genomics. Unlock the power of scale with high-throughput multiomic single cell profiling. [pages.10xgenomics.com/rs/446-PBO-704/images/10x\\_LIT000126\\_Product\\_Sheet\\_Single\\_Cell\\_GE\\_HT\\_Kit.pdf](https://pages.10xgenomics.com/rs/446-PBO-704/images/10x_LIT000126_Product_Sheet_Single_Cell_GE_HT_Kit.pdf). Accessed January 3, 2022.
9. Illumina. QC and rebalancing of single-cell gene expression libraries using the iSeq 100 System. [illumina.com/content/dam/illumina-marketing/documents/systems/iseq/single-cell-library-qc-app-note-770-2019-029.PDF](https://illumina.com/content/dam/illumina-marketing/documents/systems/iseq/single-cell-library-qc-app-note-770-2019-029.PDF). Accessed January 3, 2022.
10. Illumina. NextSeq 1000 and NextSeq 2000 single-cell RNA sequencing solution. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/nextseq-1000-2000-single-cell-app-note-m-gl-00478/nextseq-1000-2000-single-cell-app-note-m-gl-00478.pdf](https://illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/nextseq-1000-2000-single-cell-app-note-m-gl-00478/nextseq-1000-2000-single-cell-app-note-m-gl-00478.pdf). Accessed February 1, 2022.
11. 10x Genomics. Sequencing Metrics & Base Composition of Single Cell 3' v3.1 Dual Index Libraries. [support.10xgenomics.com/single-cell-gene-expression/sequencing/doc/technical-note-sequencing-metrics-and-base-composition-of-single-cell-3-v31-dual-index-libraries](https://support.10xgenomics.com/single-cell-gene-expression/sequencing/doc/technical-note-sequencing-metrics-and-base-composition-of-single-cell-3-v31-dual-index-libraries). Accessed January 3, 2022.
12. Illumina. BaseSpace Sequence Hub Demo Data. [basespace.illumina.com/datacentral](https://basespace.illumina.com/datacentral). Accessed January 4, 2022.
13. 10x Genomics. Interpreting Cell Ranger Web Summary Files for Single Cell Gene Expression Assays. [support.10xgenomics.com/single-cell-gene-expression/sequencing/doc/technical-note-interpreting-cell-ranger-web-summary-files-for-single-cell-gene-expression-assays](https://support.10xgenomics.com/single-cell-gene-expression/sequencing/doc/technical-note-interpreting-cell-ranger-web-summary-files-for-single-cell-gene-expression-assays). Accessed January 3, 2022.

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製品	カタログ番号
NovaSeq 6000 SP Reagent Kit v1.5 (100 cycles)	20028401
NovaSeq 6000 S1 Reagent Kit v1.5 (100 cycles)	20028319
NovaSeq 6000 S2 Reagent Kit v1.5 (100 cycles)	20028316
NovaSeq 6000 S4 Reagent Kit v1.5 (200 cycles)	20028313
NextSeq 1000/2000 P2 Reagents v3 (100 cycles)	20046811
NextSeq 2000 P3 Reagents (100 cycles)	20040559
NextSeq 500/550 Mid Output Kit v2.5 (150 cycles)	20024904
NextSeq 500/550 High Output Kit v2.5 (150 cycles)	20024907

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© 2023 Illumina, Inc. All rights reserved.