

InfiniumTM CytoSNP-850K v1.3 BeadChip

Comprehensive coverage of
cytogenomic-relevant genes
for constitutional and cancer
research applications

- Designed based on input from the international cytogenomics community and peer-reviewed literature
- Includes ~ 850K SNPs with 15× redundancy and enriched coverage for 3262 dosage-sensitive genes
- Offers high signal-to-noise ratios using long 50-mer SNP probes for high target specificity
- Delivers reproducible data with a broad range of sample types, including FFPE samples

illumina®

Introduction

Structural and numerical genetic variations are known to influence the presentation of a phenotype. Accurate profiling of chromosomal aberrations, such as duplications, deletions, unbalanced rearrangements, and copy-neutral absence of heterozygosity (AOH) events, is crucial for studies associated with cancer and genetic disease. The Infinium CytoSNP-850K v1.3 BeadChip (Figure 1) uses proven Infinium assay chemistry to offer the sensitivity and broad coverage needed to understand the impact of these genetic variants.

The Infinium CytoSNP-850K v1.3 BeadChip includes input on cytogenetic-relevant genes from the international community for constitutional and cancer research applications. The gene list contains updated content from the International Collaboration for Clinical Genomics (ICCG)¹ and the Cancer Cytogenomics Microarray Consortium (CCMC),² providing a comprehensive view of cytogenomic activity.

The Infinium CytoSNP-850K v1.3 BeadChip provides robust performance across a broad range of sample types, including formalin-fixed, paraffin-embedded (FFPE) tissue samples. The processed arrays can be scanned using the iScan™ System, HiScan™ System, or NextSeq™ 550 System with high reproducibility (Table 1). For exceptional data analysis, CytoSNP- 850K BeadChips are compatible with the BlueFuse™ Multi Software.

Inclusive coverage

The Infinium CytoSNP-850K v1.3 BeadChip contains ~ 850K empirically selected single nucleotide polymorphisms (SNPs) spanning across the genome. This high density of SNPs enables high-resolution analysis for discovery of meaningful chromosomal aberrations (Table 2). Intelligent design, based on the latest updates from peer-reviewed literature, provides enriched coverage for 3262 genes of known relevance in both constitutional and cancer research applications (Table 3).

For greater sensitivity, the Infinium CytoSNP-850K v1.3 BeadChip takes advantage of proven Infinium assay technology. The use of 50-mer SNP probes promotes high specificity to the target sequence. This facilitates enhanced identification of low-level mosaics³ and




Figure 1: Infinium CytoSNP-850K v1.3 BeadChip—The 8-sample Infinium CytoSNP-850K v1.3 BeadChip provides cytogenomic-relevant genes for constitutional and cancer studies with content from ICCG and CCMC.

accurate breakpoint estimation for both copy number variations (CNVs) and AOH. High 15× bead redundancy increases the signal-to-noise ratio (SNR), facilitating identification of CNV and AOH calls with high confidence. Accurate CNV calls can be made with as few as 10 consecutive probes, demonstrating the high performance level of the Infinium assay.

Powerful analysis software

BlueFuse Multi Software uses optimized, validated algorithms to detect copy number change and AOH automatically, enabling rapid, accurate molecular cytogenetics analysis. The software offers rich genomic annotation, a centralized database of historical case findings, and comprehensive reporting of identified variants to facilitate data interpretation. Developed in coordination with the clinical genomics community, the intuitive BlueFuse Multi Software interface provides an integrated framework for analyzing data for molecular cytogenetic applications.

Table 1: Infinium CytoSNP-850K BeadChip product information

| Feature | Description | | |
|-----------------------------|-------------------|--------------------|-----------------------|
| Species | Human | | |
| No. of samples per BeadChip | 8 | | |
| DNA input requirement | 200 ng | | |
| Assay chemistry | Infinium HD Super | | |
| SNP replicates | 15x | | |
| No. of SNPs to call CNV | 10 | | |
| Instrument support | iScan System | NextSeq 550 System | |
| Total no. of markers | 848,511 | 848,511 | |
| Sample throughput per week | 960 | 128 | |
| Scan time per sample | 5 min | 5 min | |
| Data performance | iScan System | NextSeq 550 System | Product specification |
| Call rate | 99.87% | 99.91% | > 98% |
| Reproducibility | 99.99% | 99.99% | >99% |
| Log R deviation | 0.1005 | 0.11 | < 0.20 |

Table 2: Infinium CytoSNP-850K v1.3 BeadChip coverage

| Probe spacing | Average distance |
|--------------------------------|------------------|
| Targeted region probe spacing | ~ 1 kb |
| Targeted region resolution | ~ 10 kb |
| Genomic backbone probe spacing | ~ 5 kb |
| Overall probe spacing | ~ 1.8 kb |
| Overall effective resolution | ~ 18 kb |

Table 3: Infinium CytoSNP-850K v1.3 BeadChip marker information^a

| Marker categories ^b | iScan System |
|--------------------------------|--------------|
| Total no. of markers | 848,511 |
| RefSeq genes | 451,304 |
| RefSeq +/- 10 kb | 522,349 |
| ADME genes | 14,147 |
| ADME +/- 10 kb | 17,419 |
| COSMIC genes | 419,775 |
| HLA markers | 5126 |
| HLA genes | 294 |
| GO genes | 139,256 |
| Exonic regions | 70,851 |
| Promoter regions | 24,891 |
| X chromosome markers | 29,837 |
| Y chromosome markers | 1229 |
| PAR/homologous markers | 984 |

a. Values are obtained from the Assay manifest file. Variations are due to different manifests/product files required to process the BeadChip for each instrument.
 b. Compared against the human genome issue hg-19 reference genome.
 Abbreviations: ADME, adsorption, distribution, metabolism, excretion; COSMIC, catalog of somatic mutations in cancer; MHC, major histocompatibility complex; HLA, human leukocyte antigen; PAR, pseudoautosomal region

Learn more

Infinium CytoSNP-850K v1.3 BeadChip, [illumina.com/techniques/popular-applications/cytogenomics.html](https://www.illumina.com/techniques/popular-applications/cytogenomics.html)

BlueFuse Multi Software, [illumina.com/clinical/clinical_informatics/bluefuse.html](https://www.illumina.com/clinical/clinical_informatics/bluefuse.html)

Ordering information

| Product | Catalog no. |
|--|-------------|
| Infinium CytoSNP-850K v1.3 BeadChip Kit (8 samples) | 20025643 |
| Infinium CytoSNP-850K v1.3 BeadChip Kit (16 samples) | 20025644 |
| Infinium CytoSNP-850K v1.3 BeadChip Kit (48 samples) | 20025645 |
| Infinium CytoSNP-850K v1.3 BeadChip Kit (96 samples) | 20025646 |

References

1. Clinical Genome Resource. www.clinicalgenome.org. Published July 2018. Accessed December 12, 2022.
2. Cancer Genomics Consortium. www.cancergenomics.org. Published July 2018. Accessed December 12, 2022.
3. Conlin LK, Thiel BD, Bonnemann CG, et al. [Mechanisms of mosaicism, chimerism and uniparental disomy identified by single nucleotide polymorphism array analysis](#). *Hum Mol Genet.* 2010;19(7):1263–1275. doi: 10.1093/hmg/ddq003.

illumina®

1.800.809.4566 toll-free (US) | +1.858.202.4566 tel
techsupport@illumina.com | www.illumina.com

© 2023 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.
M-GL-01507 v1.0